





Tera Harz

TC-80DP

TC-80DP is a high abrasion resistance and strength. It is easy to finish and polish and can be stained with all types of composite staining kits. It enables to produce a long-term prosthesis with very high strength.

- High Flexural Strengh 220 Mpa ~ 360 Mpa
- High Tensile and Impact Strength.
- Nontoxic, Biocompatible.
- For dental(C&B) and laminate.











Comparison to general C&B

Name	Tera Harz (TC-80DP)	Varseo Smile Crown plus	C&B MFH	CrownTec
	(Graphy)	(BEGO)	(NextDent)	(Saremco)
Color	A1,A2,A3 (The rest use colorants)	A1,A2,A3	A1,A2,A3	Permanent Restorations A1,A2,A3
Viscosity	2,200 cPS	$2,500 \sim 6,000 \text{ cPS}$	-	-
Bi-axial Flexural Strength	≥ 360 MPa			
Flexural strength	220 MPa	116 Mpa	107 MPa	130 MPa
Flexural modulus	4,500 MPa	4,090 MPa	-	4,000 MPa
Elongation	<10 %	<4 %	<4 %	<4 %
Wavelength	405 nm	405 nm	405 nm	405 nm

Clinical

All-on-6













Clinical Cases





Clinical Cases ediate implant placement & Loading












Clinical

Cases









Clinical







